

SPEX[®] | CHIRON

PER- AND POLYFLUOROALKYL SUBSTANCE (PFAS) STANDARDS

ENVIRONMENTAL

Explore Our Innovative & Expanding Portfolio for PFAS Testing

PFAS testing can be highly challenging which is why we make it easier with our convenient portfolio of PFAS reference materials. Advance environmental and health safety with standards you can trust.

Product Features:

- Over 390 standards for your unique applications

PFAS RMs:

- Various singles available
- Solutions in assorted solvents
- Useful for various matrices
- Custom standards of other PFAS compounds available on demand



[CHIRON.NO/EN](https://chiron.no/en)

+47 73 87 44 90

SALES@CHIRON.NO

ZeptoMetrix[®]
an antylia scientific company

Olefines

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15665.7-50-ME	Perfluorohept-1-ene	Perfluoro-1-heptene	355-63-5	50 µg/mL	methanol

Telomeric olefines

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15585.8-50-ME	1H,1H,2H-Perfluoro-1-octene	(Perfluorohexyl) ethene	25291-17-2	50 µg/mL	methanol
15586.10-50-ME	1H,1H,2H-Perfluoro-1-decene	(Perfluorooctyl) ethene	21652-58-4	50 µg/mL	methanol
15755.12-50-ME	1H,1H,2H-Perfluoro-1-dodecene	(Perfluorodecyl) ethylene	30389-25-4	50 µg/mL	methanol
15678.14-50-ME	1H,1H,2H-Perfluoro-1-tetradecene	(Perfluorododecyl) ethylene	67103-05-3	50 µg/mL	methanol
15668.16-50-ME	1H,1H,2H-Perfluoro-1-hexadecene	Perfluorotetradecyl ethylene	104564-28-5	50 µg/mL	methanol

Carboxylic acids and salts

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15118.2-50-ME	Perfluoroethanoic acid	PFEA; TFA	76-05-1	50 µg/mL	MeOH with 4 mole eq. NaOH
15141.2-100-ME	Trifluoroacetate-1-13C, sodium salt	TFA-13C1	286425-32-9	100 µg/mL	methanol
2190.3-50-ME	Perfluoropropionic acid, sodium salt	PFPrA-Na; Pentafluoropropionic acid sodium salt	378-77-8	50 µg salt/mL	methanol
2810.4-50-ME	Perfluorobutyric acid	PFBA; Heptafluorobutyric acid	375-22-4	50 µg/mL	MeOH with 4 mole eq. NaOH
2189.4-50-ME	Perfluorobutyric acid, sodium salt	PFBA-Na; Heptafluorobutyric acid sodium salt	2218-54-4	50 µg conj. base/mL	methanol
2819.5-50-ME	<i>n</i> -Perfluoropentanoic acid	PFPA; <i>n</i> -Nonafluoropentanoic acid	2706-90-3	50 µg/mL	MeOH with 4 mole eq. NaOH
2590.6-50-ME	Perfluorohexanoic acid	PFHxA	307-24-4	50 µg/mL	MeOH with 4 mole eq. NaOH
2821.7-50-ME	<i>n</i> -Perfluoroheptanoic acid	PFHpA	375-85-9	50 µg/mL	MeOH with 4 mole eq. NaOH
2820.7-50-ME	Perfluorocyclohexanecarboxylic acid	PFCyHxA	374-88-9	50 µg/mL	MeOH with 4 mole eq. NaOH
2042.8-10-AN	<i>n</i> -Perfluorooctanoic acid	PFOA	335-67-1	10 µg/mL	acetonitrile
2042.8-100-AN	<i>n</i> -Perfluorooctanoic acid	PFOA	335-67-1	100 µg/mL	acetonitrile
2042.8-10-ME	<i>n</i> -Perfluorooctanoic acid	PFOA	335-67-1	10 µg/mL	MeOH with 4 mole eq. NaOH
2042.8-50-ME	<i>n</i> -Perfluorooctanoic acid	PFOA	335-67-1	50 µg/mL	MeOH with 4 mole eq. NaOH
2822.8-50-ME	<i>n</i> -Perfluorooctanoic acid, ammonium salt	PFOA-NH ₄	3825-26-1	50 µg salt/mL	methanol
2191.8-50-ME	Perfluorooctanoic acid, sodium salt	PFOA-Na	335-95-5	50 µg salt/mL	methanol
2715.9-50-ME	Perfluorononanoic acid	PFNA	375-95-1	50 µg/mL	MeOH with 4 mole eq. NaOH
2832.9-50-ME	Perfluoro-3,5,5-trimethylhexanoic acid	PF3,5,5TriMeHxA	238403-51-5	50 µg/mL	MeOH with 4 mole eq. NaOH
2823.10-50-ME	<i>n</i> -Perfluorodecanoic acid	PFDA	335-76-2	50 µg/mL	MeOH with 4 mole eq. NaOH
2824.10-50-ME	Perfluoro-3,7-dimethyloctanoic acid	PF3,7DiMeOA	172155-07-6	50 µg/mL	MeOH with 4 mole eq. NaOH
2874.11-50-ME	Perfluoroundecanoic acid, <i>n</i> -isomer (major)	PFUndA	2058-94-8	50 µg/mL	MeOH with 4 mole eq. NaOH
2826.12-50-ME	<i>n</i> -Perfluorododecanoic acid	PFDodA	307-55-1	50 µg/mL	MeOH with 4 mole eq. NaOH
11367.13-50-ME	Perfluorotridecanoic acid, <i>n</i> -isomer (major)	PFTrDA; PFTriDA	72629-94-8	50 µg/mL	MeOH with 4 mole eq. NaOH
16021.13-50-ME	<i>n</i> -Perfluorotridecanoic acid	PFTrDA; PFTriDA	72629-94-8	50 µg/mL	MeOH with 4 mole eq. NaOH
2827.14-50-ME	Perfluorotetradecanoic acid, <i>n</i> -isomer (major)	PFTedA	376-06-7	50 µg/mL	MeOH with 4 mole eq. NaOH
15262.15-50-ME	Perfluoropentadecanoic acid (techn.)	PFPeDA	141074-63-7	50 µg/mL	MeOH with 4 mole eq. NaOH
2828.16-50-ME	Perfluorohexadecanoic acid	PFHxdA	67905-19-5	50 µg/mL	MeOH with 4 mole eq. NaOH
2829.18-50-ME	Perfluorooctadecanoic acid	PFOdA	16517-11-6	50 µg/mL	MeOH with 4 mole eq. NaOH

Dicarboxylic acids

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15379.9-50-ME	Perfluorononanedioic acid	Tetradecafluoroazelaic acid	23453-64-7	50 µg/mL	MeOH with 4 mole eq. NaOH
15380.12-50-ME	Perfluorododecanedioic acid	Eicosafluorododecanedioic acid	865-85-0	50 µg/mL	MeOH with 4 mole eq. NaOH

Carboxylic esters

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15684.11-50-ME	Dimethyl perfluoroazelate	-	22116-90-1	50 µg/mL	methanol
15682.8-50-AN	Butyl perfluorobutanoate	Butyl perfluorobutyrate	1559-07-5	50 µg/mL	acetonitrile
15689.8-50-ME	Methyl perfluoroheptanoate	-	14312-89-1	50 µg/mL	methanol
12744.9-50-ME	Methyl perfluorooctanoate	-	376-27-2	50 µg/mL	methanol
12745.10-50-AN	Ethyl perfluorooctanoate	-	3108-24-5	50 µg/mL	acetonitrile
15688.11-50-ME	Methyl perfluorodecanoate	-	307-79-9	50 µg/mL	methanol
15686.14-50-AN	Ethyl perfluorododecanoate	Tricosafluorododecanoic acid ethyl ester	158607-41-1	50 µg/mL	acetonitrile

Telomeric carboxylic esters

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15574.8-50-ME	Ethyl 2,3,3,4,4,5,5,6,6,7,7,7-dodecafluoroheptanoate	-	2338-77-4	50 µg/mL	methanol

Telomeric acids, FTCAs

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15034.8-50-ME	2H,2H-Perfluorooctanoic acid	6:2 FTCA; H2PFOA	53826-12-3	50 µg/mL	MeOH with 4 mole eq. NaOH
8548.10-50-ME	2H,2H-Perfluorodecanoic acid	8:2 FTCA; H2PFDA	27854-31-5	50 µg/mL	MeOH with 4 mole eq. NaOH
8547.11-50-ME	2H,2H,3H,3H-Perfluoroundecanoic acid	8:3 FTCA; H4PFUnA; H-PFUDa	34598-33-9	50 µg/mL	MeOH with 4 mole eq. NaOH
15619.12-50-ME	2H,2H,3H,3H-Perfluorododecanoic acid	10:2 FTCA; FDEA	53826-13-4	50 µg/mL	MeOH with 4 mole eq. NaOH
2830.6-50-ME	2H,2H,3H,3H-Perfluorohexanoic acid	3:3 FTCA	356-02-5	50 µg/mL	MeOH with 4 mole eq. NaOH
14847.8-50-ME	2H,2H,3H,3H-Perfluorooctanoic acid	5:3 FTCA	914637-49-3	50 µg/mL	MeOH with 4 mole eq. NaOH
2831.9-50-ME	2H,2H,3H,3H-Perfluorononanoic acid	6:3 FTCA	27854-30-4	50 µg/mL	MeOH with 4 mole eq. NaOH
15417.10-50-ME	2H,2H,3H,3H-Perfluorodecanoic acid	7:3 FTCA; FHpPA	812-70-4	50 µg/mL	MeOH with 4 mole eq. NaOH

Unsaturated telomeric acid, FTUCAs and unsaturated acids

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14801.6-50-ME	4,5,5,5-Tetrafluoro-4-(trifluoromethyl)-2-pentenoic acid	FTMeUPA	243139-64-2	50 µg/mL	MeOH with 4 mole eq. NaOH
14802.6-50-ME	4-(Trifluoromethylhexafluoro)-2-pentenoic acid	PFMeUPA	239795-58-5	50 µg/mL	MeOH with 4 mole eq. NaOH
2833.6-50-ME	2H,3H-Perfluoro-2-hexenoic acid	3:3 FTUCA; 3:3 fluorotelomer unsaturated carboxylic acid	37759-76-5	50 µg/mL	MeOH with 4 mole eq. NaOH
15264.8-50-ME	2H-Perfluoro-2-octenoic acid	6:2 FTUCA; 6:2 fluorotelomer unsaturated carboxylic acid	70887-88-6	50 µg/mL	MeOH with 4 mole eq. NaOH
14173.10-50-ME	2H-Perfluoro-2-decenoic acid	8:2 FTUCA; 8:2 fluorotelomer unsaturated carboxylic acid	70887-84-2	50 µg/mL	MeOH with 4 mole eq. NaOH
14474.12-50-ME	2H-Perfluoro-2-dodecenoic acid	10:2 FTUCA; 10:2 fluorotelomer unsaturated carboxylic acid	70887-94-4	50 µg/mL	MeOH with 4 mole eq. NaOH

XH perfluorocarboxylic acids

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15050.2-50-ME	2H-Perfluoroethanoic acid	2H PFEA; HPETBA; 2H-Difluoroethanoic acid	381-73-7	50 µg/mL	MeOH with 4 mole eq. NaOH
15051.3-50-ME	3H-Perfluoropropanoic acid	3H PFPrA; HPFPrA; 3H-Tetrafluoropropanoic acid	756-09-2	50 µg/mL	MeOH with 4 mole eq. NaOH
15052.4-50-ME	4H-Perfluorobutanoic acid	4H PFBA; HPFBA; 4H-Hexafluorobutanoic acid	679-12-9	50 µg/mL	MeOH with 4 mole eq. NaOH
2834.5-50-ME	5H-Perfluoropentanoic acid	5H PFPeA; HPFPeA; 5H-Octafluoropentanoic acid	376-72-7	50 µg/mL	MeOH with 4 mole eq. NaOH
2835.7-50-ME	7H-Perfluoroheptanoic acid	7H PFHPA; HPFHpA	1546-95-8	50 µg/mL	MeOH with 4 mole eq. NaOH
15053.8-50-ME	8H-Perfluorooctanoic acid	8H PFOA; HPFOA	13973-14-3	50 µg/mL	MeOH with 4 mole eq. NaOH
2836.9-50-ME	9H-Perfluorononanoic acid	9H PFNA; HPFNA	76-21-1	50 µg/mL	MeOH with 4 mole eq. NaOH
2825.11-50-ME	11H-Perfluoroundecanoic acid	11H PFUndA; HPFUndA	1765-48-6	50 µg/mL	MeOH with 4 mole eq. NaOH

XH perfluorosulfonates

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15744.4-50-ME	1,1,2,2,3,3,4,4-Octafluorobutane-1-sulphonic acid	-	70259-86-8	50 µg/mL	methanol

Telomeric alcohols

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14416.2-50-ME	2,2,2-Trifluoroethanol	1:1 FTOH	75-89-8	50 µg/mL	methanol
14829.3-50-ME	2,2-Difluoropropan-1-ol	3H,3H,3H 2:1 FTOH	33420-52-9	50 µg/mL	methanol
14828.3-50-ME	3,3,3-Trifluoropropan-1-ol	1:2 FTOH	2240-88-2	50 µg/mL	methanol
14827.3-50-ME	2,2,3,3,3-Pentafluoropropan-1-ol	2:1 FTOH	422-05-9	50 µg/mL	methanol
14831.4-50-ME	1H,1H,2H,2H,3H,3H-Perfluorobutan-1-ol	1:3 FTOH	461-18-7	50 µg/mL	methanol
2839.4-50-ME	1H,1H,2H,2H-Perfluorobutan-1-ol	2:2 FTOH	54949-74-5	50 µg/mL	methanol
14830.4-50-ME	1H,1H-Perfluorobutan-1-ol	3:1 FTOH	375-01-9	50 µg/mL	methanol
8802.6-50-ME	1H,1H,2H,2H-Perfluorohexan-1-ol	4:2 FTOH	2043-47-2	50 µg/mL	methanol
2841.6-50-ME	1H,1H,2H,2H,3H,3H-Perfluorohexan-1-ol	3:3 FTOH	679-02-7	50 µg/mL	methanol
2842.6-50-ME	1H,1H-Perfluorohexan-1-ol	5:1 FTOH	423-46-1	50 µg/mL	methanol
15708.7-50-ME	1H,1H,2H,2H-Perfluoroheptan-1-ol	5:2s FTOH	185689-57-0	50 µg/mL	methanol
2843.7-50-ME	1H,1H-Perfluoroheptan-1-ol	6:1 FTOH	375-82-6	50 µg/mL	methanol
2147.8-50-ME	1H,1H,2H,2H-Perfluorooctan-1-ol	6:2 FTOH	647-42-7	50 µg/mL	methanol
2145.8-10-ME	1H,1H-Perfluorooctan-1-ol	7:1 FTOH	307-30-2	10 µg/mL	methanol
2145.8-50-ME	1H,1H-Perfluorooctan-1-ol	7:1 FTOH	307-30-2	50 µg/mL	methanol
15741.9-50-ME*	1H,1H,2H,2H-Perfluorooctan-1-ol	7:2s FTOH	755-02-2	50 µg/mL	methanol
2844.9-50-ME	1H,1H,2H,2H,3H,3H-Perfluorononan-1-ol	6:3 FTOH	80806-68-4	50 µg/mL	methanol
2847.9-50-ME	1H,1H,2H,2H-Perfluoro-7-methyloctan-1-ol	7Me 6:1 FTOH	20015-46-7	50 µg/mL	methanol
2846.9-50-ME	1H,1H,9H-Perfluorononan-1-ol	9H 8:1 FTOH	376-18-1	50 µg/mL	methanol
2845.9-50-ME	1H,1H-Perfluorononan-1-ol	8:1 FTOH	423-56-3	50 µg/mL	methanol
15549.9-50-ME	1H,1H,1H,2H-Perfluorononan-2-ol	7:2 FTOH	24015-83-6	50 µg/mL	methanol
2851.10-50-ME	1H,1H,10H,10H-Perfluoro-1,10-decanediol	1:8:1 FTdiOH	754-96-1	50 µg/mL	methanol
2849.10-50-ME	1H,1H,2H,2H-Perfluorodecan-1-ol	8:2 FTOH	678-39-7	50 µg/mL	methanol
2848.10-50-ME	1H,1H-Perfluorodecan-1-ol	9:1 FTOH	307-37-9	50 µg/mL	methanol
2850.10-50-ME	1H,1H-Perfluoro-3,7-dimethyloctan-1-ol	3,7DiMe 7:1 FTOH	232587-50-7	50 µg/mL	methanol
2854.11-50-ME	1H,1H,2H,2H-Perfluoro-9-methyldecan-1-ol	9Me 8:2 FTOH	31200-98-3	50 µg/mL	methanol
2853.11-50-ME	1H,1H,11H-Perfluoroundecan-1-ol	11H 10:1 FTOH	307-70-0	50 µg/mL	methanol
2852.11-50-ME	1H,1H-Perfluoroundecan-1-ol	10:1 FTOH	307-46-0	50 µg/mL	methanol
2767.12-50-ME	1H,1H,2H,2H-Perfluorododecan-1-ol	10:2 FTOH	865-86-1	50 µg/mL	methanol
15557.14-50-IO	1H,1H,2H,2H-Perfluorododecyl acetate	10:2 FTOAc	37858-05-2	50 µg/mL	isooctane
3821.13-50-ME	1H,1H,13H-Perfluorotridecan-1-ol	13H 12:1 FTOH	423-72-3	50 µg/mL	methanol
12220.14-50-ME	1H,1H,2H,2H-Perfluorotetradecan-1-ol	12:2 FTOH	39239-77-5	50 µg/mL	methanol
2855.14-50-ME	1H,1H-Perfluorotetradecan-1-ol	13:1 FTOH	15622-57-8	50 µg/mL	methanol
2856.16-50-ME	1H,1H-Perfluorohexadecan-1-ol	15:1 FTOH	216144-94-4	50 µg/mL	methanol
2857.18-50-ME	1H,1H-Perfluorooctadecan-1-ol	17:1 FTOH	242142-82-1	50 µg/mL	methanol

* Minimum order of 5 required

Halogenides (iodides)

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15677.3-50-AN	Perfluoropropyl iodide	Heptafluoropropyl iodide; 66409; R21711	754-34-7	50 µg/mL	acetonitrile
15670.4-50-AN	Perfluorobutyl iodide	Nonafluoro- <i>n</i> -butyl iodide	423-39-2	50 µg/mL	acetonitrile
15676.5-50-AN	Perfluoropentyl iodide	Perfluoroamyl iodide	638-79-9	50 µg/mL	acetonitrile
15583.6-50-ME	Perfluorohexyl iodide	Tridecafluorohexyl iodide	355-43-1	50 µg/mL	methanol
15671.7-50-AN	Perfluoro- <i>n</i> -heptyl iodide	Perfluoro-1-iodoheptane	335-58-0	50 µg/mL	acetonitrile
15673.8-50-AN	Perfluorooctyl iodide	Heptadecafluoro-1-iodooctane; Telomer AN	507-63-1	50 µg/mL	acetonitrile
15582.12-50-ME	Perfluorododecyl iodide	Pentacosafuorododecyl iodide	307-60-8	50 µg/mL	methanol

Telomeric halogenides (iodides)

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15669.6-50-AN	1H,1H,2H,2H-Perfluorohexyl iodide	4:2 Fluorotelomer iodide; 4:2 FTI; Perfluorobutylethyl iodide	2043-55-2	50 µg/mL	acetonitrile
15584.8-50-ME	1H,1H,2H,2H-Perfluorooctyl iodide	6:2 Fluorotelomer iodide; 6:2 FTI; Perfluorohexylethyl iodide	2043-57-4	50 µg/mL	methanol
15666.10-50-AN	1H,1H,2H,2H-Perfluorodecyl iodide	8:2 Fluorotelomer iodide; 8:2 FTI; Perfluorooctylethyl iodide	2043-53-0	50 µg/mL	acetonitrile
15480.12-50-ME	1H,1H,2H,2H-Perfluorododecyl iodide	10:2 Fluorotelomer iodide; 10:2 FTI; Perfluorodecylethyl iodide	2043-54-1	50 µg/mL	methanol
10114.14-50-ME	1H,1H,2H,2H-Perfluorotetradecyl iodide	12:2 Fluorotelomer iodide; 12:2 FTI; Perfluorolaurylethyl iodide	30046-31-2	50 µg/mL	methanol

Ethers

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15681.14-50-ME	2H-Perfluoro-5,8,11-trimethyl-3,6,9,12-tetraoxapentadecane	Du Pont E 4; Freon O4; E 4 (oxygenator)	26738-51-2	50 µg/mL	methanol
15692.6-50-ME	1H,1H,2'H-Perfluorodipropyl ether		1000-28-8	50 µg/mL	methanol
14513.5-50-ME	Trifluoro(heptafluoro-1-propoxy) ethylene	PPVE; Cheminox FPVE; Perfluoro(propylvinylether)	1623-05-8	50 µg/mL	methanol
14514.6-50-ME	1,1,2,2,3,3-Hexafluoro-1-(trifluoromethoxy)-3-[(1,2,2-trifluorovinyl)oxy]propane	MV 31	40573-09-9	50 µg/mL	methanol
14398.9-K-IO	2-Trifluoromethyl-3-ethoxydodecafluorohexane	Novac HFE 7500; HFE 7500	297730-93-9	1000 µg/mL	isooctane

Telomeric alcohol internal standards

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
11202.6-50-ME	1D,1D,2H,2H-1-Hydroxyperfluorohexane	4:2 FTOH-d2	N/A	50 µg/mL	methanol
11203.8-50-ME	1D,1D,2H,2H-1-Hydroxyperfluorooctane	6:2 FTOH-d2	647-42-7 (unlabelled)	50 µg/mL	methanol
9825.10-50-ME	1D,1D,2H,2H-1-Hydroxyperfluorodecane	8:2 FTOH-d2	N/A	50 µg/mL	methanol
10089.12-50-ME	1D,1D,2H,2H-1-Hydroxyperfluorododecane	10:2 FTOH-d2	N/A	50 µg/mL	methanol

Allylic alcohols telomers

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
2858.6-50-ME	1H,1H,2H,3H-Perfluorohex-2-en-1-ol	Allylic 3:3 FTOH	37759-88-9	50 µg/mL	methanol
2859.9-50-ME	1H,1H,2H,3H-Perfluoronon-2-en-1-ol (cis/trans isomer mix)	Allylic 6:3 FTOH	83707-01-1/83706-97-2 (cis/trans)	50 µg/mL	methanol

Esters of telomeric alcohols including internal standards

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
10698.11-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -octyl acrylate	6:2 FTAC; 6:2 Perfluorotelomer acrylate; 6:2 FTO acrylate	17527-29-6	50 µg/mL	isooctane
10698.11-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -octyl acrylate	6:2 FTAC; 6:2 Perfluorotelomer acrylate; 6:2 FTO acrylate	17527-29-6	50 µg/mL	acetonitrile
11204.11-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -octyl acrylate-d3	6:2 FTAC-d3; 6:2 FTO acrylic acid-d3	17527-29-6 (unlabelled)	50 µg/mL	isooctane
11204.11-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -octyl acrylate-d3	6:2 FTAC-d3; 6:2 FTO acrylic acid-d3	17527-29-6 (unlabelled)	50 µg/mL	acetonitrile
9228.13-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -decyl acrylate	8:2 FTAC; 8:2 FTO acrylic acid; Repellion 120 C1	27905-45-9	50 µg/mL	isooctane
9228.13-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -decyl acrylate	8:2 FTAC; 8:2 FTO acrylic acid; Repellion 120 C1	27905-45-9	50 µg/mL	acetonitrile
9520.13-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -decyl acrylate-d3	8:2 FTAC-d3; 8:2 FTO acrylic acid-d3	N/A	50 µg/mL	isooctane
9520.13-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -decyl acrylate-d3	8:2 FTAC-d3; 8:2 FTO acrylic acid-d3	N/A	50 µg/mL	acetonitrile
10699.15-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -dodecyl acrylate	10:2 FTAC; 10:2 Fluorotelomer acrylate	17741-60-5	50 µg/mL	isooctane
10699.15-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -dodecyl acrylate	10:2 FTAC; 10:2 Fluorotelomer acrylate	17741-60-5	50 µg/mL	acetonitrile

Esters of telomeric alcohols including internal standards

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
11205.15-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -dodecyl acrylate-d3	10:2 FTAC-d3; 10:2 FTO acrylic acid-d3	N/A	50 µg/mL	isooctane
11205.15-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -dodecyl acrylate-d3	10:2 FTAC-d3; 10:2 FTO acrylic acid-d3	N/A	50 µg/mL	acetonitrile
11206.12-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -octyl methacrylate	6:2 FTMAC; 6:2 FTO methacrylic acid	2144-53-8	50 µg/mL	isooctane
11206.12-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -octyl methacrylate	6:2 FTMAC; 6:2 FTO methacrylic acid	2144-53-8	50 µg/mL	acetonitrile
11207.12-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -octyl methacrylate-d5	6:2 FTMAC-d5; 6:2 FTO methacrylic acid-d5	2144-53-8 (unlabelled)	50 µg/mL	isooctane
11207.12-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -octyl methacrylate-d5	6:2 FTMAC-d5; 6:2 FTO methacrylic acid-d5	2144-53-8 (unlabelled)	50 µg/mL	acetonitrile
15687.13-50-AN	1H,1H,9H-Perfluorononyl methacrylate	1H,1H,9H-Hexadecafluorononyl methacrylate	1841-46-9	50 µg/mL	acetonitrile
9519.14-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -decyl methacrylate	8:2 FTMAC; 8:2 FTO methacrylic acid	1996-88-9	50 µg/mL	isooctane
9519.14-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -decyl methacrylate	8:2 FTMAC; 8:2 FTO methacrylic acid	1996-88-9	50 µg/mL	acetonitrile
9521.14-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -decyl methacrylate-d5	8:2 FTMAC-d5; 8:2 FTO methacrylic acid-d5	N/A	50 µg/mL	isooctane
9521.14-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -decyl methacrylate-d5	8:2 FTMAC-d5; 8:2 FTO methacrylic acid-d5	N/A	50 µg/mL	acetonitrile
11208.16-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -dodecyl methacrylate	10:2 FTMAC; 10:2 FTO methacrylic acid	2144-54-9	50 µg/mL	isooctane
11208.16-50-AN	1H,1H,2H,2H-Perfluoro- <i>n</i> -dodecyl methacrylate	10:2 FTMAC; 10:2 FTO methacrylic acid	2144-54-9	50 µg/mL	acetonitrile
11209.16-50-IO	1H,1H,2H,2H-Perfluoro- <i>n</i> -dodecyl methacrylate-d5	10:2 FTMAC-d5; 10:2 FTO methacrylic acid-d5	N/A	50 µg/mL	isooctane

Aldehydes

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14898.7-50-ME	1H,7H-Perfluoroheptane-1,1-diol	1H,7H-Perfluoro-1-heptaldehyde hydrate	812-87-3	50 µg/mL	methanol

Alkylamines

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15694.6-50-ME	1H,1H-Perfluorohexylamine		355-34-0	50 µg/mL	methanol

Alkylamides

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15119.2-50-ME	Perfluoroethanamide	Trifluoroacetamide	354-38-1	50 µg/mL	methanol
15120.3-50-ME	Perfluoropropanamide	Pentafluoropropanamide	354-76-7	50 µg/mL	methanol
2837.4-50-ME	Perfluorobutyramide	Heptafluorobutyramide	662-50-0	50 µg/mL	methanol
15138.6-50-ME	Perfluorohexanamide	Undecafluorohexanamide	335-54-6	50 µg/mL	methanol
2838.8-50-ME	Perfluorooctanamide	Pentadecafluorooctanamide	423-54-1	50 µg/mL	methanol
15184.8-50-ME	Perfluorooctanamidine	Pentadecafluorooctanimidamide	307-31-3	50 µg/mL	methanol
15770.9-50-ME	<i>N</i> -[3-(Dimethylamino)propyl]perfluoro- <i>n</i> -butanamide	-	64790-29-0	50 µg/mL	methanol
15771.10-50-ME	<i>N</i> -[3-(Dimethylamino)propyl]perfluoro- <i>n</i> -pentanamide	-	153339-13-0	50 µg/mL	methanol

Sulfonic acids

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15140.2-50-ME	Perfluoroethanesulfonic acid, potassium salt	PFES-K	2837-92-5	50 µg conj. base/mL	methanol
14849.3-50-ME	Perfluoropropanesulfonic acid, sodium salt	PFPrS-Na	359868-82-9	50 µg conj. base/mL	methanol
2041.4-10-ME	Perfluorobutanesulfonic acid	PFBS	375-73-5	10 µg/mL	methanol
2041.4-50-ME	Perfluorobutanesulfonic acid	PFBS	375-73-5	50 µg/mL	methanol
2719.4-50-ME	Perfluorobutanesulfonic acid, potassium salt	PFBS-K	29420-49-3	50 µg conj. base/mL	methanol
15263.5-50-ME	Perfluoropentanesulfonic acid	PFPeS	2706-91-4	50 µg/mL	methanol
14797.5-50-ME	Perfluoropentanesulfonic acid, sodium salt	PFPeS-Na	630402-22-1	50 µg conj. base/mL	methanol
8579.6-50-ME	Perfluorohexanesulfonic acid	PFHxS	355-46-4	50 µg/mL	methanol
8581.6-50-ME	Perfluorohexanesulfonic acid, potassium salt	PFHxS-K	3871-99-6	50 µg conj. base/mL	methanol
12433.7-50-ME	Perfluoroheptanesulfonic acid	PFHpS	375-92-8	50 µg/mL	methanol
11156.7-50-ME	Perfluoroheptanesulfonic acid, potassium salt	PFHpS-K	60270-55-5	50 µg conj. base/mL	methanol

Sulfonic acids

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
12467.7-50-ME	Perfluoroheptanesulfonic acid, sodium salt	PFHpS-Na	21934-50-9	50 µg conj. base/mL	methanol
2037.8-50-ME	Perfluorooctanesulfonic acid, techn. isomer mixture	PFOS isomer mix.	1763-23-1 (for the <i>n</i> -isomer)	50 µg/mL	methanol
14862.8-50-ME	Perfluorooctanesulfonic acid, techn. isomer mixture, purified	PFOS isomer mix.	1763-23-1 (for the <i>n</i> -isomer)	50 µg/mL	methanol
14910.8-50-ME	Perfluorooctanesulfonic acid (<i>n</i> -isomer enriched)	PFOS, <i>n</i> -isomer enriched	1763-23-1	50 µg/mL	methanol
14863.8-50-ME	<i>n</i> -Perfluorooctanesulfonic acid	PFOS	1763-23-1	50 µg/mL	methanol
2193.8-50-ME	Perfluorooctanesulfonic acid, potassium salt	PFOS-K	2795-39-3	50 µg conj. base/mL	methanol
14717.30-50-ME	Didecyldimethylammonium perfluorooctane sulfonate	PFOS-N	251099-16-8	50 µg/mL	methanol
12746.9-50-ME	Methyl perfluorooctane sulfonate	-	71417-25-9	50 µg/mL	methanol
12747.10-50-ME	Ethyl perfluorooctane sulfonate	-	71417-26-0	50 µg/mL	methanol
2192.8-50-ME	Perfluoro-4-ethylcyclohexanesulfonic acid, potassium salt	PFECHS-K	335-24-0	50 µg conj. base/mL	methanol
14799.9-50-ME	Perfluorononanesulfonic acid, sodium salt	PFNS-Na	98789-57-2	50 µg conj. base/mL	methanol
11157.10-50-ME	Perfluorodecanesulfonic acid	PFDS	335-77-3	50 µg/mL	methanol
11309.10-50-ME	Perfluorodecanesulfonic acid, sodium salt	PFDS-Na	2806-15-7	50 µg conj. base/mL	methanol
14963.12-50-ME	Perfluorododecanesulfonic acid	PFD _o DS	79780-39-5	50 µg/mL	methanol
14947.12-50-MEX5	Perfluorododecanesulfonic acid, sodium salt	PFD _o DS-Na	1260224-54-1	50 µg conj. base/mL	methanol

Methyl sulfonates and disulfonates

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15961.13-K-ME	1H,1H,2H,2H-Perfluoro (methanesulfonyl)dodecane	10:2 FTSO ₂ Me; 10:2 Fluorotelomer methylsulfone	N/A	1000 µg/mL	methanol
15962.15-K-ME	1H,1H,2H,2H-Perfluoro (methanesulfonyl)tetradecane	12:2 FTSO ₂ Me; 12:2 Fluorotelomer methylsulfone	N/A	1000 µg/mL	methanol
15963.13-K-ME	1H,1H,2H,2H-Perfluoro (chloromethanesulfonyl)dodecane	10:2 FTSO ₂ MeCl; 10:2 Fluorotelomer chloromethylsulfone	N/A	1000 µg/mL	methanol
15964.15-K-ME	1H,1H,2H,2H-Perfluoro (chloromethanesulfonyl)tetradecane	12:2 FTSO ₂ MeCl; 12:2 Fluorotelomer chloromethylsulfone	N/A	1000 µg/mL	methanol
15965.16-K-ME	Bis(6:2Perfluorooctyl)sulfone	Bis(6:2 FT)SO ₂	N/A	1000 µg/mL	methanol

Telomeric sulfonic acid, FTSA

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15137.4-50-ME	1H,1H,2H,2H-Perfluorobutanesulfonic acid	2:2 FTSA; 2:2-FtS	149246-63-9	50 µg/mL	methanol
12748.6-50-ME	1H,1H,2H,2H-Perfluorohexanesulfonic acid	4:2 FTSA; 4:2-FtS	757124-72-4	50 µg/mL	methanol
2559.8-50-ME	1H,1H,2H,2H-Perfluorooctanesulfonic acid	6:2 FTSA; 6:2-FtS	27619-97-2	50 µg/mL	methanol
12749.10-50-ME	1H,1H,2H,2H-Perfluorodecanesulfonic acid	8:2 FTSA; 8:2-FtS	39108-34-4	50 µg/mL	methanol
12750.12-50-ME	1H,1H,2H,2H-Perfluorododecanesulfonic acid	10:2 FTSA; 10:2-FtS	120226-60-0	50 µg/mL	methanol

Sulfonamides

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15117.2-50-ME	Perfluoroethanesulfonamide	FESA	78491-70-0	50 µg/mL	methanol
15139.3-50-ME*	Perfluoropropansulfonamide	FPrSA	152894-03-6	50 µg/mL	methanol
14170.4-50-ME	<i>n</i> -Perfluorobutanesulfonamide	FBSA	30334-69-1	50 µg/mL	methanol
14902.6-50-ME	<i>n</i> -Perfluorohexanesulfonamide	FHSA	41997-13-1	50 µg/mL	methanol
14903.6-50-ME	Perfluorohexanesulfonamide, techn. Isomer mixture	FHSA	41997-13-1	50 µg/mL	methanol
2043.8-10-AN	Perfluorooctanesulfonamide, techn., <i>n</i> -isomer (major)	FOSA; PFOSA	754-91-6	10 µg/mL	acetonitrile
2043.8-50-AN	Perfluorooctanesulfonamide, techn., <i>n</i> -isomer (major)	FOSA; PFOSA	754-91-6	50 µg/mL	acetonitrile
2043.8-100-ME	Perfluorooctanesulfonamide, techn., <i>n</i> -isomer (major)	FOSA; PFOSA	754-91-6	100 µg/mL	methanol

* Minimum order of 5 required

Sulfonamides, *N*-methyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14171.5-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorobutanesulfonamide	MeFBSA	68298-12-4	50 µg/mL	methanol
14952.5-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorobutanesulfonamide-d3	MeFBSA-d3	N/A	50 µg/mL	methanol
14904.7-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorohexanesulfonamide	<i>N</i> -MeFHSA, <i>n</i> -isomer	68259-15-4	50 µg/mL	methanol
15037.7-50-ME	<i>N</i> -Methylperfluorohexanesulfonamide, isomer mixture	<i>N</i> -MeFHSA, isomer mixture	68259-15-4 as <i>n</i> -isomer	50 µg/mL	methanol
14905.7-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorohexanesulfonamide-d3	<i>N</i> -MeFHSA-d3, <i>n</i> -isomer	N/A	50 µg/mL	methanol
11932.9-50-AN	<i>N</i> -Methyl- <i>n</i> -perfluorooctanesulfonamide	<i>N</i> -MeFOSA, <i>n</i> -isomer	31506-32-8	50 µg/mL	acetonitrile
11932.9-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorooctanesulfonamide	<i>N</i> -MeFOSA, <i>n</i> -isomer	31506-32-8	50 µg/mL	methanol
13580.9-50-AN	<i>N</i> -Methylperfluorooctanesulfonamide, isomer mix	<i>N</i> -MeFOSA, isomer mixture	N/A	50 µg/mL	acetonitrile
13260.9-50-AN	<i>N</i> -Methyl- <i>n</i> -perfluorooctanesulfonamide-d3	<i>N</i> -MeFOSA-d3, <i>n</i> -isomer	936109-37-4	50 µg/mL	acetonitrile
13579.9-50-AN	<i>N</i> -Methylperfluorooctanesulfonamide-d3, isomer mix	<i>N</i> -MeFOSA-d3, isomer mix	N/A	50 µg/mL	acetonitrile
13586.10-50-AN	<i>N,N</i> -Dimethyl- <i>n</i> -perfluorooctanesulfonamide	<i>N,N</i> -Me2FOSA, <i>n</i> -isomer	213181-78-3	50 µg/mL	acetonitrile

Sulfonamides, *N*-ethyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14953.6-50-AN	<i>N</i> -Ethyl- <i>n</i> -perfluorobutanesulfonamide	<i>N</i> -EtFBSA	40630-67-9	50 µg/mL	acetonitrile
14954.6-50-AN	<i>N</i> -Ethyl- <i>n</i> -perfluorobutanesulfonamide-d5	<i>N</i> -EtFBSA-d5	N/A	50 µg/mL	acetonitrile
15809.6-50-ME	<i>N</i> -(2-Hydroxyethyl)perfluoro- <i>n</i> -butanesulfonamide	FBSAE	34454-99-4	50 µg/mL	methanol
14955.8-50-AN	<i>N</i> -Ethyl- <i>n</i> -perfluorohexanesulfonamide	<i>N</i> -EtFHSA	87988-56-5	50 µg/mL	acetonitrile
14956.8-50-AN	<i>N</i> -Ethyl- <i>n</i> -perfluorohexanesulfonamide-d5	<i>N</i> -EtFHSA-d5	N/A	50 µg/mL	acetonitrile
11933.10-50-AN	<i>N</i> -Ethyl- <i>n</i> -perfluorooctanesulfonamide	<i>N</i> -EtFOSA, <i>n</i> -isomer; Sulfluramid, <i>n</i> -isomer	4151-50-2	50 µg/mL	acetonitrile
13587.10-50-AN	<i>N</i> -Ethylperfluorooctanesulfonamide, isomer mix	<i>N</i> -EtFOSA, isomer mix; Sulfluramid, isomer mix	N/A	50 µg/mL	acetonitrile
13412.10-50-AN	<i>N</i> -Ethyl- <i>n</i> -perfluorooctanesulfonamide-d5	<i>N</i> -EtFOSA-d5, <i>n</i> -isomer; Sulfluramid-d5, <i>n</i> -isomer	936109-40-9	50 µg/mL	acetonitrile
13588.10-50-AN	<i>N</i> -Ethylperfluorooctanesulfonamide-d5, isomer mix	<i>N</i> -EtFOSA-d5, isomer mix; Sulfluramid-d5, isomer mix	936109-40-9 (for <i>n</i> -isomer)	50 µg/mL	acetonitrile

Sulfonamides, *N*-propyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15551.11-50-AN	<i>N</i> -Propyl- <i>n</i> -perfluorooctanesulfonamide	<i>N</i> -PrFOSA, <i>n</i> -isomer	2266-83-3	50 µg/mL	acetonitrile
15699.11-50-AN	<i>N</i> -Propyl- <i>n</i> -perfluorooctanesulfonamide-d7	<i>N</i> -PrFOSA-d7, <i>n</i> -isomer	N/A	50 µg/mL	acetonitrile

Sulfonamides, *N*-methyl, *N*-hydroxyethyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15253.7-50-AN	<i>N</i> -(2-Hydroxyethyl)- <i>N</i> -methylperfluorobutanesulfonamide	<i>N</i> -MeFBSE	34454-97-2	50 µg/mL	acetonitrile
15254.7-50-AN	<i>N</i> -(2-Hydroxyethyl)- <i>N</i> -methylperfluorobutanesulfonamide-d7	<i>N</i> -MeFBSE-d7	N/A	50 µg/mL	acetonitrile
15255.9-50-AN	<i>N</i> -(2-Hydroxyethyl)- <i>N</i> -methylperfluorohexanesulfonamide	<i>N</i> -MeFHxSE	68555-75-9	50 µg/mL	acetonitrile
15256.9-50-AN	<i>N</i> -(2-Hydroxyethyl)- <i>N</i> -methylperfluorohexanesulfonamide-d7	<i>N</i> -MeFHxSE-d7	N/A	50 µg/mL	acetonitrile
10979.11-50-AN	<i>N</i> -(2-Hydroxyethyl)- <i>N</i> -methylperfluorooctanesulfonamide	<i>N</i> -MeFOSE	24448-09-7	50 µg/mL	acetonitrile
10979.11-50-ME	<i>N</i> -(2-Hydroxyethyl)- <i>N</i> -methylperfluorooctanesulfonamide	<i>N</i> -MeFOSE	24448-09-7	50 µg/mL	methanol
14846.11-50-ME	<i>N</i> -(2-Hydroxyethyl)- <i>N</i> -methylperfluorooctanesulfonamide-d7	<i>N</i> -MeFOSE-d7	1265205-95-5	50 µg/mL	methanol

Sulfonamides, *N*-ethyl, *N*-hydroxyethyl

Volume: 1 mL

15257.8-50-AN	<i>N</i> -Ethyl- <i>N</i> -(2-hydroxyethyl)perfluorobutanesulphonamide	<i>N</i> -EtFBSE	34449-89-3	50 µg/mL	acetonitrile
15258.8-50-AN	<i>N</i> -Ethyl- <i>N</i> -(2-hydroxyethyl)perfluorobutanesulphonamide-d9	<i>N</i> -EtFBSE-d9	N/A	50 µg/mL	acetonitrile
15259.10-50-AN	<i>N</i> -Ethyl- <i>N</i> -(2-hydroxyethyl)perfluorohexanesulphonamide	<i>N</i> -EtFHxSE	34455-03-3	50 µg/mL	acetonitrile
15352.10-50-AN	<i>N</i> -Ethyl- <i>N</i> -(2-hydroxyethyl)perfluorohexanesulphonamide-d9	<i>N</i> -EtFHxSE-d9	N/A	50 µg/mL	acetonitrile
11931.12-50-AN	<i>N</i> -Ethyl- <i>N</i> -(2-hydroxyethyl)perfluorooctanesulphonamide	<i>N</i> -EtFOSE	1691-99-2	50 µg/mL	acetonitrile
11931.12-50-ME	<i>N</i> -Ethyl- <i>N</i> -(2-hydroxyethyl)perfluorooctanesulphonamide	-	1691-99-2	50 µg/mL	methanol
14848.12-50-ME	<i>N</i> -Ethyl- <i>N</i> -(2-hydroxyethyl)perfluorooctanesulphonamide-d9	<i>N</i> -EtFOSE-d9	1691-99-2 (unlabelled)	50 µg/mL	methanol

Sulfonamides, *N*-propyl, *N*-hydroxyethyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15552.13-50-AN	<i>N</i> -Propyl- <i>N</i> -(2-hydroxyethyl) perfluorooctanesulphonamide	<i>N</i> -PrFOSE	N/A	50 µg/mL	acetonitrile
15698.13-50-AN	<i>N</i> -Propyl- <i>N</i> -(2-hydroxyethyl) perfluorooctanesulphonamide-d11	<i>N</i> -PrFOSE-d11	N/A	50 µg/mL	acetonitrile

Sulfonylglycine

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14901.6-50-ME	<i>N</i> -Perfluorobutanesulfonylglycine	FBSAA	347872-22-4	50 µg/mL	methanol
14957.8-50-ME	<i>N</i> -Perfluorohexanesulfonylglycine	FHxSAA	1003193-99-4	50 µg/mL	methanol
15778.8-50-ME	Perfluorohexanesulfonylglycine, isomer mixture	FHxSAA, isomer mix	1003193-99-4 (<i>n</i> -isomer)	50 µg/mL	methanol
14526.10-50-AN	<i>N</i> -Perfluorooctanesulfonylglycine	FOSAA	2806-24-8	50 µg/mL	acetonitrile
15772.12-50-ME	<i>N</i> -(Carboxymethyl)- <i>N</i> -[(heptadecafluorooctyl)sulfonyl]glycine	-	251359-47-4	50 µg/mL	methanol
15773.10-50-ME	<i>N</i> -(Carboxymethyl)- <i>N</i> -[(tridecafluorohexyl)sulfonyl]glycine	-	N/A	50 µg/mL	methanol

Sulfonylglycine, *N*-methyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14172.7-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorobutanesulfonylglycine	<i>N</i> -MeFBSAA	159381-10-9	50 µg/mL	methanol
14949.7-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorobutanesulfonylglycine-d3	<i>N</i> -MeFBSAA-d3	N/A	50 µg/mL	methanol
14950.9-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorohexanesulfonylglycine	<i>N</i> -MeFHxSAA	715646-50-7	50 µg/mL	methanol
14951.9-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorohexanesulfonylglycine-d3	<i>N</i> -MeFHxSAA-d3	N/A	50 µg/mL	methanol
13257.11-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorooctanesulfonylglycine	<i>N</i> -MeFOSAA, <i>n</i> -isomer	2355-31-9	50 µg/mL	methanol
13258.11-50-ME	<i>N</i> -Methyl- <i>n</i> -perfluorooctanesulfonylglycine-d3	<i>N</i> -MeFOSAA-d3, <i>n</i> -isomer	1400690-70-1	50 µg/mL	methanol
13585.11-50-ME	<i>N</i> -Methylperfluorooctanesulfonylglycine-d3, isomer mix	<i>N</i> -MeFOSAA-d3, isomer mix	N/A	50 µg/mL	methanol

Sulfonylglycine, *N*-ethyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14958.8-50-ME	<i>N</i> -Ethyl- <i>n</i> -perfluorobutanesulfonylglycine	<i>N</i> -EtFBSAA	68957-33-5	50 µg/mL	methanol
14959.8-50-ME	<i>N</i> -Ethyl- <i>n</i> -perfluorobutanesulfonylglycine-d5	<i>N</i> -EtFBSAA-d5	N/A	50 µg/mL	methanol
14960.10-50-ME	<i>N</i> -Ethyl- <i>n</i> -perfluorohexanesulfonylglycine	<i>N</i> -EtFHSAA	68957-32-4	50 µg/mL	methanol
14961.10-50-ME	<i>N</i> -Ethyl- <i>n</i> -perfluorohexanesulfonylglycine-d5	<i>N</i> -EtFHSAA-d5	N/A	50 µg/mL	methanol
13259.12-50-ME	<i>N</i> -Ethyl- <i>n</i> -perfluorooctanesulfonylglycine	<i>N</i> -EtFOSAA, <i>n</i> -isomer	2991-50-6	50 µg/mL	methanol
13590.12-50-ME	<i>N</i> -Ethylperfluorooctanesulfonylglycine, isomer mix	<i>N</i> -EtFOSAA, isomer mix	N/A	50 µg/mL	methanol
13413.12-50-ME	<i>N</i> -Ethyl- <i>n</i> -perfluorooctanesulfonylglycine-d5	<i>N</i> -EtFOSAA-d5, <i>n</i> -isomer	1265205-97-7	50 µg/mL	methanol
13591.12-50-ME	<i>N</i> -Ethylperfluorooctanesulfonylglycine-d5, isomer mix	<i>N</i> -EtFOSAA-d5, isomer mix	N/A	50 µg/mL	methanol

Sulfonylglycine, *N*-carboxymethyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15717.8-50-ME	<i>N</i> -(Carboxymethyl)- <i>N</i> -[(1,1,2,2,3,3,4,4,4-nonafluorobutyl)sulfonyl]glycine	FBSEE-DA	1268835-43-3	50 µg/mL	methanol

Other sulfonamides/-amines

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15807.14-50-ME	3-[[3-(Dimethylamino)propyl][(tridecafluorohexyl)sulfonyl]amino]-1-propanesulfonic acid	SPrAmPr-FHxSA; FH-SA- <i>N</i> -3-dimethylaminopropyl- <i>N</i> -3-propanesulfonic acid	38850-60-1	50 µg/mL	methanol
15808.9-50-ME	3-[[[(Perfluorohexyl)sulfonyl]amino]propanesulfonic acid	SPr-FHxSA; FH-SA- <i>N</i> -3-propanesulfonic acid	2089108-64-3	50 µg/mL	methanol
15707.15-50-ME	6:2 Fluorotelomermercaptoalkylamido sulfonate	6:2 FTSAS	N/A	50 µg/mL	methanol

Thioethers and sulfones

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15805.11-50-ME	3-[(1H,1H,2H,2H-Perfluoro- <i>n</i> -octyl)thio]propanoic acid	6:2 FTThPrA; 6:2 Perfluorohexylthio 3-propanoic acid	149339-57-1	50 µg/mL	methanol
15806.11-50-ME	3-[(1H,1H,2H,2H-Perfluoro- <i>n</i> -octyl)sulfonyl]propanoic acid	6:2 FTSO2PrA; 6:2 Perfluorohexylsulfone 3-propanoic acid	2089110-04-1	50 µg/mL	methanol

Amphoteric sulfonamides and amides

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15799.10-50-ME	3-[[Perfluoro- <i>n</i> -butyl)sulfonyl]amino]- <i>N,N,N</i> -trimethyl-1-propanaminium bromide	PFBSAmS; FB-SA- <i>N</i> -(<i>N,N,N</i> -trimethyl-3-propanaminium) bromide	N/A	50 µg conj. base/mL	methanol
15800.12-50-ME	3-[[Perfluoro- <i>n</i> -hexyl)sulfonyl]amino]- <i>N,N,N</i> -trimethyl-1-propanaminium bromide	PFHxSAmS; FH-SA- <i>N</i> -(<i>N,N,N</i> -trimethyl-3-propanaminium) bromide	N/A	50 µg conj. base/mL	methanol
15801.14-50-ME	3-[[Perfluoro- <i>n</i> -octyl)sulfonyl]amino]- <i>N,N,N</i> -trimethyl-1-propanaminium bromide	PFOSAmS; FO-SA- <i>N</i> -(<i>N,N,N</i> -trimethyl-3-propanaminium) bromide	73149-44-7	50 µg conj. base/mL	methanol
15745.11-50-ME	<i>N</i> -[3-(Dimethylamino)propyl]perfluoro- <i>n</i> -hexanamide	PFHxAAM	153339-14-1	50 µg/mL	methanol
15802.9-50-ME	<i>N</i> -[3-(Dimethylamino)propyl]perfluoro- <i>n</i> -butanesulfonamide	AmPr-FBSA; FB-SA- <i>N</i> -(<i>N,N,N</i> -trimethyl-3-propanamine)	68555-77-1	50 µg/mL	methanol
15803.11-50-ME	<i>N</i> -[3-(Dimethylamino)propyl]perfluoro- <i>n</i> -hexanesulfonamide	AmPr-FHxSA; FH-SA- <i>N</i> -(<i>N,N,N</i> -trimethyl-3-propanamine)	50598-28-2	50 µg/mL	methanol
15804.13-50-ME	<i>N</i> -[3-(Dimethylamino)propyl]perfluoro- <i>n</i> -octylsulfonamide	AmPr-FOSA; FO-SA- <i>N</i> -(<i>N,N,N</i> -trimethyl-3-propanamine)	13417-01-1	50 µg/mL	methanol
15947.13-50-ME	Capstone product A	DPOSA	80475-32-7	50 µg/mL	methanol
15774.17-50-ME	8:2 Fluorotelomer thioether amido sulfonic acid	8:2 FtTAoS	755698-73-8	50 µg/mL	methanol
15775.13-50-ME	4:2 Fluorotelomer thioether amido sulfonic acid	4:2 FtTAoS	1333933-57-5	50 µg/mL	methanol
15776.15-50-ME	6:2 Fluorotelomer sulfonyl amido sulfonic acid	6:2 FtSO2AoS	1911606-13-7	50 µg/mL	methanol

Sulfonyl acrylate, *N*-methyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15332.10-50-IO	<i>N</i> -Methylperfluorobutanesulphonamidoethyl acrylate	MeFBSEA	67584-55-8	50 µg/mL	isooctane
15333.10-50-IO	<i>N</i> -Methylperfluorobutanesulphonamidoethyl acrylate-d7	MeFBSEA-d7	N/A	50 µg/mL	isooctane
11857.12-50-IO	<i>N</i> -Methylperfluorohexanesulphonamidoethyl acrylate	MeFHxSEA	67584-57-0	50 µg/mL	isooctane
15334.12-50-IO	<i>N</i> -Methylperfluorohexanesulphonamidoethyl acrylate-d7	MeFHxSEA-d7	N/A	50 µg/mL	isooctane
15067.14-50-IO	<i>N</i> -Methylperfluorooctanesulfonamidoethyl acrylate	MeFOSEA	25268-77-3	50 µg/mL	isooctane
15335.14-50-IO	<i>N</i> -Methylperfluorooctanesulfonamidoethyl acrylate-d7	MeFOSEA-d7	N/A	50 µg/mL	isooctane

Sulfonyl acrylate, *N*-ethyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15336.11-50-IO	<i>N</i> -Ethylperfluorobutanesulphonamidoethyl acrylate	EtFBSEA	17329-79-2	50 µg/mL	isooctane
15337.11-50-IO	<i>N</i> -Ethylperfluorobutanesulphonamidoethyl acrylate-d9	EtFBSEA-d9	N/A	50 µg/mL	isooctane
15339.13-50-IO	<i>N</i> -Ethylperfluorohexanesulphonamidoethyl acrylate	EtFHxSEA	1893-52-3	50 µg/mL	isooctane
15340.13-50-IO	<i>N</i> -Ethylperfluorohexanesulphonamidoethyl acrylate-d9	EtFHxSEA-d9	N/A	50 µg/mL	isooctane
15341.15-50-IO	<i>N</i> -Ethylperfluorooctanesulfonamidoethyl acrylate	EtFOSEA	423-82-5	50 µg/mL	isooctane
15342.15-50-IO	<i>N</i> -Ethylperfluorooctanesulfonamidoethyl acrylate-d9	EtFOSEA-d9	N/A	50 µg/mL	isooctane

Sulfonyl acrylate, *N*-propyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15553.16-50-IO	<i>N</i> -Propylperfluorooctanesulfonamidoethyl acrylate	PeFOSEA	2357-60-0	50 µg/mL	isooctane
15697.16-50-IO	<i>N</i> -Propylperfluorooctanesulfonamidoethyl acrylate-d11	PeFOSEA-d11	N/A	50 µg/mL	isooctane

Sulfonyl methacrylate, *N*-methyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15343.11-50-IO	<i>N</i> -Methylperfluorobutanesulphonamidoethyl methacrylate	MeFBSEMA	67584-59-2	50 µg/mL	isooctane
15344.11-50-IO	<i>N</i> -Methylperfluorobutanesulphonamidoethyl methacrylate-d7	MeFBSEMA-d7	N/A	50 µg/mL	isooctane
15345.13-50-IO	<i>N</i> -Methylperfluorohexanesulphonamidoethyl methacrylate	MeFHxSEMA	67584-61-6	50 µg/mL	isooctane
15346.13-50-IO	<i>N</i> -Methylperfluorohexanesulphonamidoethyl methacrylate-d7	MeFHxSEMA-d7	N/A	50 µg/mL	isooctane
15347.15-50-IO	<i>N</i> -Methylperfluorooctanesulphonamidoethyl methacrylate	MeFOSEMA	14650-24-9	50 µg/mL	isooctane
15348.15-50-IO	<i>N</i> -Methylperfluorooctanesulphonamidoethyl methacrylate-d7	MeFOSEMA-d7	N/A	50 µg/mL	isooctane

Sulfonyl methacrylate, *N*-ethyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15349.12-50-IO	<i>N</i> -Ethylperfluorobutanesulphonamidoethyl methacrylate	EtFBSEMA	67939-33-7	50 µg/mL	isooctane
15350.12-50-IO	<i>N</i> -Ethylperfluorobutanesulphonamidoethyl methacrylate-d9	EtFBSEMA-d9	N/A	50 µg/mL	isooctane
15351.14-50-IO	<i>N</i> -Ethylperfluorohexanesulphonamidoethyl methacrylate	EtFHxSEMA	67906-70-1	50 µg/mL	isooctane
15353.14-50-IO	<i>N</i> -Ethylperfluorohexanesulphonamidoethyl methacrylate-d9	EtFHxSEMA-d9	N/A	50 µg/mL	isooctane
15354.16-50-IO	<i>N</i> -Ethylperfluorooctanesulphonamidoethyl methacrylate	EtFOSEMA	376-14-7	50 µg/mL	isooctane
15355.16-50-IO	<i>N</i> -Ethylperfluorooctanesulphonamidoethyl methacrylate-d9	EtFOSEMA-d9	N/A	50 µg/mL	isooctane

Sulfonyl methacrylate, *N*-propyl

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15554.17-50-IO	<i>N</i> -Propylperfluorooctanesulphonamidoethyl methacrylate	PrFOSEMA	13285-40-0	50 µg/mL	isooctane
15696.17-50-IO	<i>N</i> -Propylperfluorooctanesulphonamidoethyl methacrylate-d11	PrFOSEMA-d11	N/A	50 µg/mL	isooctane

GenX and related

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14110.6-50-ME	2,3,3,3-Tetrafluoro-2-(heptafluoropropoxy)propionic acid	HFPO-DA; "GENX" precursor; FRD-903; PFPrOPrA	13252-13-6	50 µg/mL	methanol
14110.6-100-ME	2,3,3,3-Tetrafluoro-2-(heptafluoropropoxy)propionic acid	HFPO-DA; "GENX" precursor; FRD-903; PFPrOPrA	13252-13-6	100 µg/mL	methanol
14112.6-50-ME	2,3,3,3-Tetrafluoro-2-(heptafluoropropoxy)propionic acid ammonium salt	GenX; FRD-902; Perfluoro-2-propoxypropanoic acid ammonium salt	62037-80-3	50 µg conj. base/mL	methanol
15128.9-50-ME	Hexafluoropropylene oxide trimer acid	HFPO-TA; Perfluoro-2,5-dimethyl-3,6-dioxanonanoic acid	13252-14-7	50 µg/mL	methanol
14115.7-50-ME	Methyl 2,3,3,3-tetrafluoro-2-(heptafluoropropoxy)propionate	HFPO-DA methyl ester; GENX methyl ester	13140-34-6	50 µg/mL	methanol
14515.5-50-IO	1,1,1,2,2,3,3-Heptafluoro-3-(1,2,2,2-tetrafluoroethoxy)propane	Freon E1; Dupont E1	3330-15-2	50 µg/mL	isooctane

ADONA and related:

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14175.7-50-ME	Ammonium 4,8-dioxa-3H-perfluor-nonanoate	ADONA (ammonium salt)	958445-44-8	50 µg conj. base/mL	methanol
14176.7-50-ME	Sodium dodecafluoro-3H-4,8-dioxanonanoate	NaDONA; ADONA sodium salt	2250081-67-3	50 µg conj. base/mL	methanol
14460.7-50-ME	Dodecafluoro-3H-4,8-dioxanonanoic acid	DONA	919005-14-4	50 µg/mL	methanol
14460.7-100-ME	Dodecafluoro-3H-4,8-dioxanonanoic acid	DONA	919005-14-4	100 µg/mL	methanol

9CI-PF3ONS and 11CI-PF3OUNS and related

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
14174.8-50-ME	9-Chlorohexadecafluoro-3-oxanone-1-sulfonic acid potassium salt	9CI-PF3ONS-K; C8CI PFESA-K; F-53B-K; 6:2 Cl-PFESA-K	73606-19-6	50 µg conj. base/mL	methanol
14462.10-50-ME	11CI-PF3OUdS potassium salt	11CI-PF3OUdS-K	83329-89-9	50 µg conj. base/mL	methanol
15882.4-50-ME	Perfluoroethoxysulfonic acid, sodium salt	NVHOS-Na	1132933-86-8	50 µg conj. base/mL	methanol
15883.4-50-ME	3-(Methoxy)tetrafluoropropionic acid	MTP	93449-21-9	50 µg/mL	methanol

Perfluoroalkyl ether sulfonic acids (PFEEAS)

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15884.7-50-ME	Perfluoro-3,6-dioxa-4-methyl-7-octenesulfonic acid	PS acid; PFESA BP1; Nafion by-product 1; Nafion BP1; NBP1	29311-67-9	50 µg/mL	methanol
15885.7-50-ME*	7H-Perfluoro-4-methyl-3,6-dioxatanesulfonic acid	Hydro-PS acid; PFESA BP2; Nafion by-product 2; Nafion BP2; NBP2	749836-20-2	50 µg/mL	methanol
14722.4-50-ME	Perfluoro(2-ethoxyethane)sulfonic acid	PFEEASA	113507-82-7	50 µg/mL	methanol
15556.7-50-ME	Perfluoro(4-methyl-3,6-dioxaoct-7-ene) sodium sulfonate	PSVNa	N/A	50 µg/mL	methanol

* Minimum order of 5 required

Perfluoroalkyl ether carboxylic acids (PFEECAs)

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15129.3-50-ME	Perfluoro-2-methoxyacetic acid	PFMOAA	674-13-5	50 µg/mL	MeOH with 4 mole eq. NaOH
14721.4-50-ME	Perfluoro-3-methoxypropanoic acid	PFMPA	377-73-1	50 µg/mL	MeOH with 4 mole eq. NaOH
14056.5-50-ME	Perfluoro-4-methoxybutanoic acid	PFMOBA; PFMBA	863090-89-5	50 µg/mL	MeOH with 4 mole eq. NaOH
15932.8-50-ME	EVE acid	Perfluoro-5-methyl-4,7-dioxanon-8-enoic acid	69087-46-3	50 µg/mL	MeOH with 4 mole eq. NaOH
15933.8-50-ME	Hydro-EVE acid	8H-Perfluoro-5-methyl-4,7-dioxanonanoic acid	773804-62-9	50 µg/mL	MeOH with 4 mole eq. NaOH
15413.5-50-ME	Perfluoro-3,6-dioxaheptanoic acid	NFDHA; 3,6-OPFHpA; Nonafluoro-3,6-dioxaheptanoic acid	151772-58-6	50 µg/mL	MeOH with 4 mole eq. NaOH
15413.5-100-ME	Perfluoro-3,6-dioxaheptanoic acid	NFDHA; 3,6-OPFHpA; Nonafluoro-3,6-dioxaheptanoic acid	151772-58-6	100 µg/mL	MeOH with 4 mole eq. NaOH
15719.6-50-ME	Perfluoro-3,6-dioxaoctanoic acid	-	80153-82-8	50 µg/mL	MeOH with 4 mole eq. NaOH
15720.6-50-ME*	Ammonium perfluoro-3,6-dioxaoctanoate	-	908020-52-0	50 µg/mL	methanol

* Please inquire

Hexanol ethers

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15704.8-50-ME	6-(2,2,2-Trifluoroethoxy)-1-hexanol	-	1341666-66-7	50 µg/mL	methanol
15705.10-50-ME	6-(2,2,3,3,4,4,4-Heptafluorobutoxy)-1-hexanol	-	1207481-10-4	50 µg/mL	methanol
15706.12-50-ME	6-[(2,2,3,3,4,4,5,5,6,6,6-undecaheptyloxy)-1-hexanol	-	N/A	50 µg/mL	methanol

Telomeric alkylsilanes

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15691.10-50-AN	1H,1H,2H,2H-Perfluorooctyldimethylchlorosilane	Chlorodimethyl[2-(perfluorohexyl)ethyl]silane	102488-47-1	50 µg/mL	acetonitrile
15675.11-50-AN	1H,1H,2H,2H-Perfluorooctyltrimethoxysilane	Tridecafluorooctyltrimethoxysilane; Dynasilan F88161; HF 81	85857-16-5	50 µg/mL	acetonitrile
15674.14-50-AN	1H,1H,2H,2H-Perfluorooctyltriethoxysilane	2-Perfluorohexylethyltriethoxysilane; Dynasilan8261; 13F (fluorosilane)	51851-37-7	50 µg/mL	acetonitrile
15667.20-50-ME	(1H,1H,2H,2H-Perfluorotetradecyl)tris(ethoxy)silane	1H,1H,2H,2H-Perfluorotetradecyltriethoxysilane	885275-56-9	50 µg/mL	methanol

Other PFAS

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15672.8-50-AN	1H,1H,2H,2H-Perfluorooctanethiol	Foralkyl EM 6; 2-(Perfluorohexyl)ethane-1-thiol	34451-26-8	50 µg/mL	acetonitrile
15683.7-50-ME	(Hydroxymethyl)perfluorocyclohexane	(Perfluorocyclohexyl)methanol	28788-68-3	50 µg/mL	methanol
15693.8-50-ME	1H,1H-Perfluoro-3,6-dioxadecan-1-ol	-	152914-73-3	50 µg/mL	methanol
15690.9-50-AN	3-Perfluorohexyl-2-iodopropan-1-ol	-	38550-44-6	50 µg/mL	acetonitrile
15716.13-50-ME	<i>N</i> -[3-(Dimethylamino)propyl]-3,3,4,4,5,5,6,6,7,7,8,8,8-tridecafluoro-1-octane-sulfonamide	M4; FS SAAM-6C	34455-22-6	50 µg/mL	methanol
15718.4-50-MEX5	1,1,2,2,3,3,4,4,4-Nonafluoro-1-butane-sulfonic acid	PFBSi	34642-43-8	50 µg/mL	methanol

PAPS

Volume: 1 mL (reconstituted)

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15280.6-50UG	4:2 MonoPAPS	Mono[2-(perfluorobutyl)ethyl] phosphate	150065-76-2	ca. 50 µg	neat
9563.8-50UG	6:2 MonoPAPS	-	57678-01-0	ca. 50 µg	neat
11340.16-50-ME	6:2 DiPAPS sodium salt	-	407582-79-0	50 µg conj. base/mL	methanol
9565.24-50UG	6:2 TriPAPS	-	165325-62-2	ca. 50 µg	neat
9391.10-50UG	8:2 MonoPAPS	-	57678-03-2	ca. 50 µg	neat
14573.10-50-ME	8:2 MonoPAPS sodium salt	-	130771-95-8 (Conjugated base)	50 µg conj. base/mL	methanol
9392.20-50UG	8:2 DiPAPS	-	678-41-1	ca. 50 µg	neat
14572.20-50-ME	8:2 DiPAPS sodium salt	-	114519-85-6	50 µg conj. base/mL	methanol
9566.12-50UG	10:2 MonoPAPS	-	57678-05-4	ca. 50 µg	neat
13264.24-50-ME	10:2 DiPAPS sodium salt	-	120433-74-1	50 µg conj. base/mL	methanol

Deuterated internal standards PAPs

Volume: 1 mL (reconstituted)

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
9928.16-50UG	6:2 DiPAPS-d4	-	N/A	ca. 50 µg	neat
9929.20-50UG	8:2 DiPAPS-d4	-	N/A	ca. 50 µg	neat
15759.24-50-ME	10:2 DiPAPS-d4 sodium salt	-	120433-74-1 (unlabelled)	50 µg conj. base/mL	methanol

FTEOs, commercial mixture

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15853.X-K-ME	Capstone® FS-30	6:2 FTEO	1640092-35-8	1000 µg/mL	methanol

FT-metacrylate

Volume: 1 mL

Part No.	Product Name	Synonym	CAS No.	Concentration	Matrix
15854.12-50-AN	Capstone® 62MA	6:2 FT Metacrylate	2144-53-8	50 µg/mL	acetonitrile

Methods Cross Reference - All 50 µg/mL in methanol, 1 mL

Carboxylic Acids

Part No.	Product Name	Synonym	CAS No.	PFAS for Key EPA Methods		
				EPA 533	EPA 537.1	EPA 1633
2810.4-50-ME	Perfluorobutyric acid	PFBA	375-22-4	x		x
2819.5-50-ME	n-Perfluoropentanoic acid	PFPeA	2706-90-3	x		x
2590.6-50-ME	Perfluorohexanoic acid	PFHxA	307-24-4	x	x	x
2821.7-50-ME	n-Perfluoroheptanoic acid	PFHpA	375-85-9	x	x	x
2042.8-50-ME	n-Perfluorooctanoic acid	PFOA	335-67-1	x	x	x
2715.9-50-ME	Perfluorononanoic acid	PFNA	375-95-1	x	x	x
2823.10-50-ME	n-Perfluorodecanoic acid	PFDA	335-76-2	x	x	x
2874.11-50-ME	Perfluoroundecanoic acid, n-isomer (major)	PFUnDA	2058-94-8	x	x	x
2826.12-50-ME	n-Perfluorododecanoic acid	PFDoDA	307-55-1	x	x	x
16021.13-50-ME	n-Perfluorotridecanoic acid	PFTTrDA; PFTiA	72629-94-8		x	x
2827.14-50-ME	Perfluorotetradecanoic acid, n-isomer (major)	PFTeDA	376-06-7		x	x

Sulfonic Acids

2041.4-50-ME	Perfluorobutanesulfonic acid	PFBS	375-73-5	x	x	x
15263.5-50-ME	Perfluoropentanesulfonic acid	PFPeS	2706-91-4	x		x
8579.6-50-ME	Perfluorohexanesulfonic acid	PFHxS	355-46-4	x	x	x
12433.7-50-ME	Perfluoroheptanesulfonic acid	PFHpS	375-92-8	x		x
14863.8-50-ME	Perfluorooctanesulfonic acid	PFOS	1763-23-1	x	x	x
14799.9-50-ME	Perfluorononanesulfonic acid	PFNS	98789-57-2			x
11157.10-50-ME	Perfluorodecanesulfonic acid	PFDS	335-77-3			x
14963.12-50-ME	Perfluorododecanesulfonic acid	PFDoS	79780-39-5			x

Telomeric sulfonic Acids

12748.6-50-ME	1H,1H,2H,2H-Perfluorohexanesulfonic acid	4:2 FTS (4:2 FTSA)	757124-72-4	x		x
2559.8-50-ME	1H,1H,2H,2H-Perfluorooctanesulfonic acid	6:2 FTS	27619-97-2	x		x
12749.10-50-ME	1H,1H,2H,2H-Perfluorodecanesulfonic acid	8:2 FTS	39108-34-4	x		x

Sulfonamides

Part No.	Product Name	Synonym	CAS No.	PFAS for Key EPA Methods		
				EPA 533	EPA 537.1	EPA 1633
16024.8-50-ME	n-Perfluorooctanesulfonamide	(P)FOSA	754-91-6	x	x	x
11932.9-50-ME	N-Methyl-n-perfluorooctanesulfonamide	NMeFOSA	31506-32-8	x		x
13259.12-50-ME	N-Ethyl-n-perfluorooctanesulfonylglycine	NEtFOSA	4151-50-2	x	x	x
10979.11-50-ME	N-(2-Hydroxyethyl)-N-methylperfluorooctane-sulfonamide	NMeFOSE	24448-09-7	x		x
11931.12-50-ME	N-Ethyl-N-(2-hydroxyethyl)perfluorooctane-sulphonamide	NEtFOSE	1691-99-2	x	x	x
13257.11-50-ME	N-Methyl-n-perfluorooctanesulfonylglycine	N-MeFOSAA	2355-31-9			x
13259.12-50-ME	N-Ethyl-n-perfluorooctanesulfonylglycine	N-EtFOSAA	2991-50-6			x

Ethers, Carboxylic acids

14175.7-50-ME	Ammonium 4,8-dioxa-3H-perfluorononanoate	ADONA	919005-14-4	(X) NH4-salt	(X) NH4-salt	(X) NH4-salt
14110.6-50-ME	2,3,3,3-Tetrafluoro-2-(heptafluoropropoxy) propionic acid	HFPO-DA	13252-13-6	x	x	x

Ethers, Carboxylic acids, other

14721.4-50-ME	Perfluoro-3-methoxypropanoic acid	PFMPA	377-73-1	x		x
15413.5-50-ME	Perfluoro-3,6-dioxaheptanoic acid	NFDHA	151772-58-6	x		x
14056.5-50-ME	Perfluoro-4-methoxybutanoic acid	PFMBA	863090-89-5	x		x

Telomeric acids

2830.6-50-ME	2H,2H,3H,3H-Perfluorohexanoic acid	3:3 FTCA	356-02-5			x
14847.8-50-ME	2H,2H,3H,3H-Perfluorooctanoic acid	5:3 FTCA	914637-49-3			x
15417.10-50-ME	2H,2H,3H,3H-Perfluorodecanoic acid	7:3 FTCA; FHpPA	812-70-4			x

Ethers, Sulfonic acids, 9Cl-11Cl etc

14462.10-50-ME	11Cl-PF3OUdS Potassium salt	11Cl-PF3OUdS-K	83329-89-9	(X) K-salt	(X) K-salt	(X) K-salt
14174.8-50-ME	9-Chlorohexadecafluoro-3-oxanone-1-sulfonic acid potassium salt	9Cl-PF3ONS K; 6:2 Cl-PFAES-K	73606-19-6	(X) K-salt	(X) K-salt	(X) K-salt

Ethers, Sulfonic acids, other

14722.4-50-ME	Perfluoro(2-ethoxyethane)sulfonic acid	PFEESA	113507-82-7	x		x
---------------	--	--------	-------------	---	--	---



CHIRON.NO/EN

+47 73 87 44 90

SALES@CHIRON.NO