

## PRODUCT DESCRIPTION:

Rubella Virus is an enveloped icosahedral Togavirus that contains a nucleocapsid and linear, non-segmented, single-stranded, positive-sense RNA.

Each frozen aliquot contains 1 mL of viral culture fluid. Rubella culture fluid is assayed by PCR (copies/mL).

## INTENDED USE:

Viral culture fluids are sold as consumable testing materials, propagation or commercialization are prohibited without prior written consent from ZeptoMetrix. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The purchase of infectious microorganisms from ZeptoMetrix requires a Material Transfer Agreement (MTA).

## BIOSAFETY:








Rubella Virus is a Biosafety Level 2 microorganism and must be used within a Biological Safety Level 2 facility or cabinet. Please consult your institution's regulations regarding the use of this organism. For a detailed discussion on biological safety see the current edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC.

## PRECAUTIONS:

- Use Universal Precautions, this product is **potentially biohazardous**.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary). Titer will be altered by multiple freeze-thaws.
- To avoid cross-contamination, use separate pipette tips for all reagents.

## RECOMMENDED STORAGE:

Viral culture fluids should be stored at -65°C or below.

	Catalog Number		Temperature Limitation
	Batch Code		Expiration Date
	For Research Use Only		Biological Risk
	Manufacturer		