

**Influenza A H1N1pdm Virus**  
**Strain: NY/02/09\*\***  
**Infectious Culture Fluid (1 mL)**

## PRODUCT DESCRIPTION:

Influenza A H1N1pdm Virus (Strain: NY/02/09) is an enveloped Orthomyxoviridae virus with a diameter of 80-120 nm. It contains a single-stranded, segmented, negative-sense RNA within a nucleocapsid.

Each frozen aliquot contains 1 mL of titered viral culture fluid.

## INTENDED USE:

Viral culture fluids are sold as consumable testing materials; propagation or commercialization is prohibited without prior written consent from ZeptoMetrix. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

**The purchase of infectious microorganisms from ZeptoMetrix requires a Material Transfer Agreement (MTA).**

**\*The purchase of live Influenza A H1N1pdm virus requires a USDA permit.**

## BIOSAFETY:

Influenza A H1N1pdm Virus is a Biosafety Level 2 microorganism and must be used within Biological Safety Level 2 facility or cabinet. Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the current edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC.

## PRECAUTIONS:








- Use Universal Precautions, this product is **potentially biohazardous**.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary). Titer will be altered by multiple freeze-thaws.
- To avoid cross-contamination, use separate pipette tips for all reagents.

## RECOMMENDED STORAGE:

Viral culture fluids should be stored at -65°C or below.

\*\*Please note that although similar in nomenclature, **this is a 2009 H1N1 pandemic Influenza strain** and does NOT correlate with the *seasonal* 2009 Influenza strains found in the FluDb.org database. For reference, the NCBI Taxon IDs for the *seasonal* Influenza strains listed in the FluDb.org database are: A/New York/01/2009 (H1N1) - 666252; B/New York/01/2009 - 664512; A/New York/02/2009 (H1N1) - 666298; and A/New York/03/2009 (H3N2) - 659637.

PI0810109CFN  
Revision: 08  
Effective Date: 09/12/2022

 REF	Catalog Number		Temperature Limitation
 LOT	Batch Code		Expiration Date
 RUO	For Research Use Only		Biological Risk
	Manufacturer		

PCA# 22-308  
Page 1 of 1