

**\*Human Metapneumovirus (hMPV-9) 9 Type A1  
Strain: IA3-2002  
Infectious Culture Fluid (1 mL)**

## PRODUCT DESCRIPTION:

Human Metapneumovirus (hMPV-9) 9 Type A1 (Strain: IA3-2002) is an enveloped virus that contains a negative-sense single-stranded RNA.

Each frozen aliquot contains 1 mL of titered viral culture fluid.

## INTENDED USE:

Viral culture fluids are sold as consumable testing materials; propagation or commercialization is prohibited without prior written consent from ZeptoMetrix. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The purchase of infectious microorganisms from ZeptoMetrix requires a Material Transfer Agreement (MTA).

## BIOSAFETY:

Human Metapneumovirus is a Biosafety Level 2 microorganism and must be used within a Biological Safety Level 2 facility or cabinet. Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the current edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC.

## PRECAUTIONS:








- Use Universal Precautions, this product is **potentially biohazardous**.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary). Titer will be altered by multiple freeze-thaws.
- To avoid cross-contamination, use separate pipette tips for all reagents.

## RECOMMENDED STORAGE:

Viral culture fluids should be stored at -65°C or below.

**\*Note: "This product is sold by ZeptoMetrix under license from ViroNovative B.V. under patent applications, including U.S. Patent Applications 10/371,099 and 10/371,122, and any patents that issue from applications related to PCT/NL02/00040 and PCT/US03/05271."**

PI0810160CF  
Revision: 06  
Effective Date: 09/19/2022

 REF	Catalog Number		Temperature Limitation
 LOT	Batch Code		Expiration Date
 RUO	For Research Use Only		Biological Risk
	Manufacturer		

www.ZeptoMetrix.com

ZeptoMetrix LLC • 878 Main Street, Buffalo, NY 14202 USA • Tel (800) 274-5487

This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

PCA# 22-308  
Page 1 of 1