

**Influenza B Virus
(Victoria Lineage)
Strain: Brisbane/46/15
Infectious Culture Fluid (1 mL)**

PRODUCT DESCRIPTION:

Influenza B Virus (Strain: Brisbane/46/15) is an enveloped virus with a diameter of 80-120 nm, and contains a single-stranded, segmented, negative-sense RNA within a nucleocapsid.

Each frozen aliquot contains 1 mL of titered viral culture fluid.

INTENDED USE:

Viral culture fluids are sold as consumable testing materials; propagation or commercialization is prohibited without prior written consent from ZeptoMetrix. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

The purchase of infectious microorganisms from ZeptoMetrix requires a Material Transfer Agreement (MTA).

BIO SAFETY:

Influenza B Virus is a Biosafety Level 2 microorganism and must be used within Biological Safety Level 2 facility or cabinet. Please consult your institution's regulations regarding the use of this organism. For a detailed discussion on biological safety see the current edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC.








PRECAUTIONS:

- Use Universal Precautions, this product is **potentially biohazardous**.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary). Titer may be altered by multiple freeze-thaws.
- To avoid cross-contamination, use separate pipette tips for all reagents.

RECOMMENDED STORAGE:

Viral culture fluids should be stored at -65°C or below.

PI0810528CF
Revision: 07
Effective Date: 02/28/2023

 REF	Catalog Number	 TEMP	Temperature Limitation
 LOT	Batch Code	 EXPI	Expiration Date
 RUO	For Research Use Only	 BIO	Biological Risk
 MFG	Manufacturer		

www.ZeptoMetrix.com

ZeptoMetrix LLC • 878 Main Street, Buffalo, NY 14202 USA • Tel (800) 274-5487

This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

PCA# 22-446
Page 1 of 1