

**SARS-Related Coronavirus 2 (SARS-CoV-2)  
Lineage JN.1; Omicron Variant<sup>1</sup>  
Isolate: USA/New York/PV96109/2023  
Culture Fluid (Heat Inactivated) (1 mL)**

Catalog Number:  
0810689CFHI

## PRODUCT DESCRIPTION:

SARS-Related Coronavirus 2 Lineage JN.1; Omicron Variant (Isolate: USA/New York/PV96109/2023) is an enveloped, positive-sense single-stranded RNA virus from the Coronaviridae family and the Betacoronaviridae genus. Severe acute respiratory syndrome-related coronavirus 2 (SARS-CoV-2), isolate hCoV-19/USA/New York/PV96109/2023 was isolated from a human nasopharyngeal swab in 2023, in New York, USA.

Under the nomenclature system introduced by GISAID (Global Initiative on Sharing All Influenza Data), SARS-CoV-2, isolate hCoV-19/USA/New York/PV96109/2023 is assigned lineage JN.1 (Pango v.4.3.1 consensus call), Omicron (BA.2-like) (Scorpio) and GISAID clade GRA using the Phylogenetic Assignment of Named Global Outbreak lineages (PANGO) tool. The complete genome of SARS-CoV-2, isolate hCoV-19/USA/New York/PV96109/2023 has been sequenced (GISAID: EPI\_ISL\_18563626). It was labeled as a variant of interest (VOI) GRA (BA.2.86+BA.2.86.\*), and first detected in Denmark/Israel/USA by the World Health Organization (WHO).

Each frozen aliquot contains 1 mL of heat inactivated viral culture fluid. The pre-inactivation titer was determined from an infectious aliquot.

Viral inactivation is verified after heat inactivation by the absence of viral growth in tissue culture-based infectivity assays.








## INTENDED USE:

Heat inactivated viral culture fluids are sold as consumable testing materials. The suitability and performance characteristics should be determined by your laboratory for each intended usage.

These products are NOT intended for use in the manufacture or processing of injectable products subject to licensure under section 351 of the Public Health Service Act or for any other product intended for administration to humans.

## FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

PI0810689CFHI  
Revision: 00  
Effective Date: 04/29/2024

|   |                       |   |                        |
|---|-----------------------|---|------------------------|
|  | Catalog Number        |  | Temperature Limitation |
|  | Batch Code            |  | Expiration Date        |
|  | For Research Use Only |  | Biological Risk        |
|  | Manufacturer          |   |                        |

www.ZeptoMetrix.com

ZeptoMetrix LLC • 878 Main Street, Buffalo, NY 14202 USA • Tel (800) 274-5487

This product was manufactured in a facility which has a Quality Management System that is ISO 13485 certified.

## BIOSAFETY:

Please consult your institution's regulations regarding the use of this product. For a detailed discussion on biological safety see the current edition of Biosafety in Microbiological and Biomedical Laboratories (BMBL), published by the CDC.

## PRECAUTIONS:

- Use Universal Precautions, this product is **potentially biohazardous**.
- Repetitive freezing and thawing is not recommended (aliquot material if necessary).
- To avoid cross-contamination, use separate pipette tips for all reagents.

## RECOMMENDED STORAGE:

Heat inactivated viral culture fluids should be stored at -65°C or below.

<sup>1</sup>The following reagent was obtained through BEI Resources, NIAID, NIH: SARS-Related Coronavirus 2, Isolate hCoV-19/USA/New York/PV96109/2023 (Lineage JN.1; Omicron Variant), NR-59693.